

At DaVinci Mobility we can make virtually any wheelchair to your exact specification — that's what we've been doing since the 1990s and we have thousands of happy customers who come back whenever they need an upgrade.

We're also experts when it comes to motor vehicle adaptations — we're

Motability accredited fitters of any motor vehicle.

We love to chat about all things mobility related to help our customers, so please get in touch to discuss your requirements.

WHAT IS BESPOKE - PAGE 4

DAVINCI LEGGERA - PAGE 6

DAVINCI XTREME - PAGE 8

DAVINCI COMPATTA - PAGE 10

DAVINCI KIDS - PAGE 12

PHOTOSHOOT - PAGE 14

CUSTOMIZE - PAGE 16

DAVINCI TRAIL RIDER - PAGE 24

DAVINCI HANDCYCLE - PAGE 28

WHAT OUR CUSTOMERS SAY - PAGE 32

CONTACT - BACKPAGE





WHAT IS BESPOKE?

Bespoke is defined as "made for a particular customer or user".

DaVinci mobility can guarantee that the chairs they sell are fully "bespoke".

Thousands of our customers had their wheelchairs made uniquely for thier size and purposes. What ever your need is, the Davinci staff will work hard to find a solution.

We are the only indpendently owned bespoke wheelchair manufactuer in the UK. This means we can create wheelchairs for you, that many of our competitors won't be able to.

In addition to this, we can also fully customize the look and feel of the wheelchair. The full spectrum of colours can be chosen for the frame and upholstery.

We pride ourselves on our customer service, we work hard to get you the chair you want. We are here to give adivce on what products will be suitable for your needs.

All Wheelchairs you see in this brochure can be adjusted to fit you. These adjustments are made after your enquiry and assessment with one of our staff.

TYPOLOGY

Seat Width
Seat Depth
Footplate Depth
Backrest Height
Backrest Angle
Camber
Front Height
Rear Height
Setting (Point of
Balance)
Front Frame Angle
Total Width
Total length
Wheelchair Weight
Maximum Load









The simple cantilever-style frame makes for a lightweight chair with very little under the seat to obstruct easy loading into a car. This design also helps to keep costs down.

Features:

- Aluminium frame (7020 Series tubing) Tough, lightweight design
- Fully adjustable custom-made upholstery
- Variable or fixed axle
- Height-adjustable backrest
- Front castors 4", 5"or 6" diameter with softroll tyres
- Large choice of powder-coated colours
- Hideaway or standard brakes
- 24", 25" or 26" rear wheels (choice of tires)





The frame shape allows for easy fitting of most popular powered and none-powered hand-crank cycles, powered front and rear wheel attachments.





The lightest rigid-framed wheelchair we make at DaVinci and one lightest wheelchairs on the market. The frame is made from a special thin-wall aluminium tubing which gives it strength yet keeps the lightness required by active everyday wheelchair users.


Features:

- Aluminium or titanium tubing
- Compact, lightweight design
- Fully-adjustable, custom-made upholstery
- Variable or fixed axle
- Height-adjustable backrest
- Front castors 4", 5", or 6"
- Large choice of colours
- Hideaway or standard brakes





Ergonomically designed to be lifted into a car by the driver and stowed either on the front or rear passenger seat. Simply remove the quick release rear wheels and you can lift the frame with ease past the steering wheel.





The way the Compatta folds keeps the rigid frame format with both rear wheels fixed on the same axle. Box frame and fixed footrest this makes the Compatta push just like a rigid frame.

Features:

- Aluminium 7020-grade tube frame
- Unique custom-built folding frame design
- Fully adjustable custom-made upholstery
- One-hand folding processHeight-adjustable backrest
- Front castors 4", 5", or 6"
- Large choice of colours
- Hideaway or standard brakes
- 24", 25" or 26"rearwheels





Once folded, the frame is very easily stowed inside a car boot or behind the passenger seat.

Ergonomically designed to be lifted into a car by the driver and stowed away easily.





The cantilever-style frame makes for a lightweight but tough frame with very little under the seat to obstruct easy loading into a car. The chair's backrest is also height adjustable with fully custom made upholstery seating.

Features:

- Aluminium Frame (7020 Series tubing) Tough,
 lightweight design
- Fully-adjustable, custom-made upholstery
- Variable axle positionHeight-adjustable backrest
- Front castors 4", 5" or 6" with soft-roll tyres
- Large choice of powder-coated colours
- Hideaway or standard brakes
- 20" 22" or 24" rear wheels (choice of tires)





The frame shape allows the fitting of most popular powered rearwheel sets. Choice of wheel size 20" 22" or 24" with hideaway or standard brakes. Additional spoke gaurds can be added with almost any custom design to give children a chance to personalise thier chair.

Si.D "GROW IT TOGETHER"

"Matsunaga" brand (SEED: seed), which is new to UK and Europe with fans of "Matsunaga" The plant "seeds" sown by using water and fertilizer to sprout, Bloom flowers and create new "seed". I want to create like this good cycle in this market with YOU.

Matsunaga's design is a simple one but very effective for optimal use, this allows you to live life the best way possible. Simple design allows YOU to grow.

Both sides from footrest area to top of frame are built without welding which makes the frame lighter and durable. Matsunaga call it "美" (Bi) FRAME. "美" means beauty.

MP original pipe, Minamo, is both rigid and light.Si:D has adopted it on the important parts for driving due to enhance strength and rigidity, which wheeling can be comfortable.

Individual purchase order form available for detailed specifica on and optons, please make sure the individual purchase order form. Specifica ons may change due to product improvement or quality improvement.



SPECIFICATION



Seat Width 300mm ∼ 420mm



Front Seat Height 10mm increments



Back Rest Height 240mm ∼ 380mm



Seat Depth 340mm ~ 420mm 10mm increments



Rear Seat Height 10mm increments



Frame Length $380 \text{mm} \sim 500 \text{mm} / 520 \text{mm} / 560 \text{mm} / 600 \text{mm}$



Camber



Maximum weight allowed



Axle Position



B-MAX AJ

Full frame adjustment for novice to veteran players

This wheelchair can be adjusted ten different ways, allowing you to try out various positions on your own. The axle tube is 40mm thick.



B-MAX

For senior players to top players, a semiadjustable frame

In B-MAX range, this model is the lightest. If you are experienced and after some adjustability, you should pick this model. The axle tube is 40mm thick. center of gravity, footrest height, wheel gap, etc. are all adjustable points.



B-MAX

This type is the one to choose if you want additional strength and rigour

The majority of the Japan national team uses this model. The axle tube is 50mm thick. Adjustable points include the footrest height and centre of gravity.
























Join The Ride





More than just a Paddle school When you paddle with us you help us to give back to the local community,

WhatSUP started to use the sport of SUP for icebreaker days for foster families, work with homeless family charities and other charities.

Now, thanks to you, we are growing fast and we can help more people get out on the water.

Collaboration is the key to success

For over 30 years the amazing team at DaVinci mobility have been making the impossible a reality and changing thousands of people's lives.

We could not be prouder to handover manufacturing and production of our disability paddle board chair to them. From us at WhatSUPuk we just want to say a huge "thank you" to the teamm, the new chair is another level.

We are always working hard behind the scenes to make paddle boarding as inclusive as we can. "Together we stand up for the community"







CUSTOMIZE YOUR CHAIR

DaVinci mobility staff work closely with thier customers to make sure the chair we create for you, fits your needs. We want you to have control of how your chair looks and behaves.

Bespoke wheelchairs made for you, that's DaVinci Mobility!







Make changes to your chair, not matter what it is we'll do our best to make it fit your needs. change frames, seats, handrims, tyres, wheels, castors, foot rests. whatever it is we'll get it working for you! Below is an example of the components of DaVinci Mobility Wheelchairs that you can configure to your liking.

SEAT + CUSHION

The seat and Cushion's Hieght, width, depth, style, embroidery and upholstery can all be made how you want.

BACKREST

The backrest can be fully customizable and taylormade to the chair users size.

HANDRIM

How ever you push your wheelchair, we'll have a handrim that suits you.



CASTORS


The castor can have many uses, Davinci offer an array of different colours and sizes.


DaVinci have a number of different frames on offer. All wheelchair frames have different attributes.

Once you've have chosen a frame, we can taylor the chair to the wheelchair user's size and the metal can be painted in any colour.

TYRE

DaVinci have a number of different tryes, what ever the purpose we'll have the tyre for you.







UPHOLSTERY





We are supported by Sunrise Medical and NHS Wheelchair Services to make sure we find the correct backrest and seating for you. Other cushions and backrests can be fitted to our wheelchairs if needed. Below are DaVinci Mobility custom made cushions and backrests. Depth, width and height for all cushion and backrests are made to your size and come in any colour you wish.


DaVinci cushions are detachable from the wheelchair seat. The cushion attaches to the wheelchair seat with velcro, this allows you to easily clean the wheelchair or change the cushion.

















2 Piece leather

Material: 2 Piece Leather

Design: Any

Colours: Any

Cost: £380

Standard
Material: Nylon
Design: Standard
Colours: Black
Cost: Free

Coloured Nylon
Material: Nylon
Design: Standard
Colours: Any
Cost: £105

<u>1 piece leather</u> Material: Nylon + Leather Design: Line or Diamond Colours: Any Cost: £330





Bergamot Moss

Thyme.

Teal.

Rosemary.

Black.

Lip.

Maderia.

Amethyst.

Imperial.

Nightshade.

Burgundy.

Totem.

Niagara.

Pewter.

Lilac Moon

Charcoal Grev.

Explore

Buffalo.

Doe.

Manilla.

Maple

Lightstone Beige

<u>Angora.</u>

Alpaca.

Birch.

Soft White.

Taupe

Rich.

Peat.

Biscuit.

Beechwood.

<u>Auburn.</u>

Highland Bark.

Cocoa.

Tulip.

Racing Yellow

Admiralty Blue

Cerestiai.

Egyptian Blue

Midnight.

Regal Blue

CUSTOMIZATION

We can fully customize the upholstery of the chair as we make our own seats and backrests. You can choose from a huge range of colours, the most popular are shown down the left hand side of the page. Below you can see other custom designs'. This includes custom logo's, we can add any logo to your seat or backrest, the possibilities are endless.





BACKREST, SEAT AND CUSHION DESIGN

Straight line Embroidery Diamond Embroidery





DaVinci backrest can have custom logos added just like the picture above. The customization ideas are endless and can add a real personal touch to your bespoke wheelchair.

loopwheels



Loopwheels Extreme

Extra Stability, Extra Speed, Extra Freedom

Loopwheels have developed a robust wheel with world-class suspension that works at higher speeds, for off-road use with a MTB tire and power attachments.

- A high-performance suspension wheel specifically designed for maximum comfort off-road when using a power attachment
- Fitted with a Mountain Bike Tire for extra grip and comfort
- Extreme comfort, extreme stability, extreme freedom: go wherever the path takes you
- Extreme ride quality and comfort on any adventure
- _

Extra Carbon, Extra Confidence, Extra Control

With 12% extra hi-tech carbon than our classic Loopwheels, Loopwheels Carbon offer a sportier ride – so you can go steadier and stronger than before.

- Designed for a sportier suspension than our Loopwheels Classic, at low speeds and primarily for manual use.
- Gorgeous carbon finish and reflective logos to help you be seen in the dark



Loopwheels Carbon



Loopwheels Multicoloured classic

Key features:

- A smoother ride on any surface
- Up to 70% vibration reduction
- With our colour options, personalise the look of your chair AND increase your comfort.
- Designed for optimum comfort on a manual wheelchair.

Key features:

- A smoother ride on any surface
- Up to 70% vibration reduction
- With our colour options, personalise the look of your chair AND increase your comfort.
- Designed for optimum comfort on a manual wheelchair, for speeds up to 7 km/h.



Loopwheels Classic





Spinergy Lite Extreme LX

Designed for people looking for a performance wheel with a unique cutting edge style of its own that offers excellent reach-through access underneath your chair. The Lite Extreme LX is aesthetically pleasing and is designed to reduce weight from your chair without sacrificing durability and performance.

Spinergy Sport Lite

Extreme SLX R-10

The SLX R-10 is a tried and true design from Spinergy with a strong and durable, stylish radial laced spoke pattern. Take the benefits of the SLX 24 spoke wheel and add the durability of a 5/8" axle. Mostly used in contact sports where the 1/2" axles tend to bend the most. These wheels deliver confidence and performance on the courts.





Spinergy Lite Extreme LX FLEXRIM

The Spinergy FLEXRIM is a revolutionary product that is changing the way people roll. With the urethane material holding the pushrim onto the rim, it allows the user to push with out ever touching the tire. The ultra-light and extremely strong material allows the pushrim to flex in and give the wheelchair user an extra inch or so when squeezing into very tight doorways. This product is unlike anything offered by any other company and a must have for everyday users of wheelchairs.

Spinergy Sport Lite Extreme SLX 18 Spoke

Get the performance you demand with Spinergy's patented 4mm PBO FIBER spokes in a wheel designed to last. This wheel is ideal for high agility sports such as Tennis and Basketball. The design with 18 spokes provides a quick maneuvering wheel that won't let you down if the playing gets rough. The Sport Lite Extreme is a great looking, stylish and performance-oriented wheel.





Spinergy Spox Everyday

Spinergy built its reputation as a market leader with the SPOX Everyday wheelchair wheel. They have been our best seller for 15 years and are built to last. Wheels are the most rewarding upgrade you can make on a wheelchair and the SPOX Everyday wheels are our entry wheel with a budget in mind.

TYRES AND CASTORS



MBL Castors Black Plastic Colour: Black



Frog Legs Castor Frog Legs Castor Frog Legs Castor Aluminium Colour: Black / Green



Aluminium Colour: Black / Blue



Aluminium Colour: Black / **Purple**



Frog Legs Castor **Aluminium** Colour: Red/ Aluminium



Frog Legs Castor **Aluminium** Colour: Black / Orange



Frog Legs Castor Aluminium Colour: Blue / Aluminium



Tire technology of the highest level is fast, reliable and agile. A perfect all-rounder for cyclists who prefer to ride with tubes.



Schwalbe Marathon Plus Evolution

The Marathon has become the absolute benchmark for puncture-proof pneumatic tyres on the wheelchair.



Schwalbe Landcruiser

The Schwalbe Landcruiser Tyre has low rolling resistance on the flats and great grip in wet conditions thanks to its raised diamond tread.



Schwalbe Crazy Bob

The complete 'FAT' beach wheel. Each wheel comes completed and ready just to fit straight onto your wheelchair.

Fitted with a Schwalbe Crazy Bob tyre which is grippy and has a wrap-around tread for added traction.



Schwalbe SpeedAir Basketball

The indoor sports specialist tyre from Schwalbe. Nimble and easy running with an asymmetrical design designed to match indoor wheelchair camber.



Schwalbe RightRun Smooth Sports

Schwalbe's entry level tire, but packed with features: Kevlar® puncture protection layer, hand-friendly 2Grip side walls, and Robust Black'n' Roll tread center compound that is non-marking neatly completes the feature list.

DaVinci wheelchairs are available in almost any colour you can imagine — from brushed titanium or polished aluminium to the full spectrum of the rainbow. You know what colour you want? We'll make it!



TRAIL RIDERS

"Power assistance at your fingertips, let the DaVinci Trailrider take the strain from pushing outdoors"

Starting From £3200













SPEC

- Available in 2 or 4 wheel drive
- 2 wheel drive has
- 2 x 1000watt brushless hub motors in the 20" rear wheels
- 4 wheels fitted with hydronic disc brakes
- 4 x 20" x 2.25" off road tyres
- 2 x 50amp motor controllers
- 2 x 48v 20amp lithium batteries
- Custom made steel tube frame powder coated in a choice of colours
- Seat fitted with Padded upholstery and cushion
 Options

OPEN A NEW DOOR

Wondering why you might need this? The DaVinci Quattro is guaranteed to deliver an unparalleled exciting and fun experience. The 4x4 has the ability to travel over all terrains making it one of the most inclusive mobility aids on the market. Unlock the door to accessibility with the DaVinci Quattro.

2 Wheel Drive - Starting From £6995 4 Wheel Drive - Starting From £9995









CICLO

"If it's cycling for leisure or to maintain cardio-vascular fitness the DaVinci 8-Speed hand-cycle is sure to meet your expectations"

Starting From £2095





Ideal for the "first-timers" when it comes to handcycling. Simply clip on the handcycle to your wheelchair, then pedal forward to get started and reverse cycle to stop. Reliability and control is provided by a shimano nexus 8 speed hub-mounted gear set for smooth gear changes and a back pedal brake for safe ans controlled stopping. With the gear shifter mounted at the top of the frame, it is very easy to select gears and keep going at the pace you want to set. If your needs change you can upgrade your handcycle from our wide range of retro fit options.

DaVinci handcycle's come in multiple options. Visit our website store or visit us instore to gather more information on which handcycle will suit your needs.

One of the lightest power assisted handcycle you can buy in the uk!

Davinci sell both manual and power assisted hancycles. The powerful yet incredibly light power assisted handcycle, comes with a LCD cycle computer and power assist selector. This completes the perfect system for when you need a little assistance.


Both the manual and power assisted handcycles can clip onto your wheelchair and you can be moving in seconds.

The Davinci handcycles are ideal for anyone who wishes to get involved in handcycling. It comes standard with 9, Sram derailleur gears and a mountain drive for propulsion. Trigger or twist grip gear shifters, disc and V brake for stopping.



- Aluminium 7020-grade tube frame
- Adjustable crank height and reach
- Fits to most rigid-frame wheelchairs (some folding frames — call to check)
- 9-speed SRAM derailleur gears with trigger or twist shifter
- 20" front wheel and choice of tire
- Powerful disk and V brakePowder-coated frame large choice of colours
- LCD cycle computer and power-assist control
- · Easy docking system with stand





SPEC

- 18speed (with mountain drive)
- 9 speed derailleur gears with thumb or twist shifter
- Mountain drive (2 speed in crank shifter)
- Front disk brake with 2 x calipers
- Quick release docking system With stand
- Fits most wheel chairs
- Choice of gear ratios
- Front wheel size 20" or 26"
- Options include rear wheel extension axle for better traction
- 10ah 36v lithium battery and charger
- All complimented with blue tooth control and display









CYCLE WITH A HELPING HAND

Are you wanting to hand-cycle but need a little support? The Davinci Veloce hand-cycle has a Keyde E-bike system which will support you as you power your way along the cycle lanes around town or the county side. Choose from 5 different levels of support as you travel around. Get active, get cycling get the Davinci Veloce!

The Battery (as shown in the picture) looks exactly like a aluminium water bottle. It's positioned down the left hand side of the hand-cycle.



The special needs bikes of Van Raam

Van Raam's special needs bikes gives people mobility. We offer a wide range of customised bicycles for people with disabilities and for people who can no longer ride a traditional bicycle or no longer have the confidence to cycle:

Individual cycle
Side-by-side tandem bikes
Tandem bikes
Wheelchair bikes
Transport bikes

Tricycles
Mobility scooter bike
Low step-through bike
Walking aid
Cycling together



Low step through bikes

The low step through bikes from Van Raam enable people to get on their bikes easier.



Walking aid

The walking aid elieves the legs by supporting the costumer's weight on the saddle.



Electric mobility scooter bike

The compact electric scooter bike is designed for people with stability problems or variable energy level.



Transport bikes

Our transport bikes are for children or adults who cannot or do not want to take part in traffic on their own or for recreational use.



Tricycles

People, who cannot ride on a regular bicycle, look for a proper and safe alternative.



Wheelchair bikes

The Van Raam wheelchair bikes are used to enable a wheelchair dependent person to travel by cycle.



Tandems

The Van Raam tandems are used to enjoy cycling together.



Side-by-side tandems

The two persons have a good view and can communicate well with each other due to the seats being next to each other.

The unique riding characteristics of every Van Raam bicycle provide a special cycling experience. With innovative and modern production techniques, and by adding options/accessories and frame adjustments, we can make a customised bike for just about anyone. The unique Van Raam pedal support is available for most of these adapted bikes.

Van Raam sells bikes all over the world, with a focus on most European countries and North America. Nothing is produced in stock; each bicycle is unique and made to order. The sale of Van Raam adapted bikes goes via the Van Raam dealers. At Van Raam in Varsseveld (the Netherlands), it is possible to make a test ride. Read more about this on the test ride-page.





'Van Raam offers sustainable and contemporary mobility to people with a disability.'





Interoduction

My name is Ellis Palmer Babe. I'm 27, a radio producer and live near Liverpool in Northwest of England. When I'm not out cycling around town, I'm supporting Manchester City, cooking pasta, or playing with my rescue cat,

So far, I'm 13,000 kilometers into my hand-cycling adventure and yet every turn of the cranks still feels liberating. Growing up, I used to look at cyclists and think "Man, I wish I wasn't disabled and I could do that". I'd see them in my cab on the way to the office and feel a pang of envy and internalized ableism.

How Ellis got into Hand-cycling

It was only two years ago when I realized I wouldn't have to miss out any longer. I was doing a coastal walk in Country Wicklow, Ireland. The battery on my wheelchair went flat "There must be a better way of getting around than this", I thought. My friend asked me if I had thought about hand-cycling and the idea kept chipping away at me.

Initially, the price tag put me off a bit and I think my parents were a bit skeptical, not just because of the cost, but also because they doubted whether I'd get much use out of it. Some 13,000 kilometer's later, it was very much worth using the overtime money I'd been saving up. It's been a lifeline!

The Hand-cycling Community

Yes, I have a great community of non-standard cyclist I can fall back on. Many of the staff at DaVinci Mobility in Liverpool who make my hand-cycles are disabled themselves - a great resource for tips and tricks around life as a wheelchair user. There's a Hand-cycling Facebook group who I can reach out to with any queries regarding, say, winter cycling, or be inspired by with different set-ups as well.























ENJOY EVERYDAY FREEDOM

How the Motability Scheme works

Exchange your qualifying mobility allowance for a brand-new car, Wheelchair Accessible Vehicle (WAVWheelchair Accessible Vehicle), scooter or powered wheelchair, so you can enjoy more freedom and independence.



Adaptations to help you drive

Adaptations can be made to many driving controls, including steering wheels, accelerators and gearsticks.



Adaptations to help you get in and out

Swivel chairs, transfer plates and hoists can make it easier for you to get in and out of your vehicle.



Adaptations to help you store equipment

You may need to take a wheelchair or something else that helps when you're out and about. Boot or roof hoists can make this easier.

POPULAR ATTACHMENTS



Jeff Gosling: Push/Pull Brake/Accelerator Hand Control

The number one choice for hand control users in the UK. Our market leading hand control is designed for drivers who are unable to operate the foot pedals. This adaptation allows a driver to operate the brake and accelerator by hand using a push/pull single lever.



Smart Lifter Boot Hoist

The Smart Lifter range of car boot hoists enable you to lift a mobility device, such as a wheelchair, scooter or powerchair, easily and quickly into your car boot at the touch of a button. This eliminates the need for the wheelchair/scooter user or family member to struggle with getting the device into a vehicle or car when wanting to venture further afield.



SmartSteer

The SmartSteer steering unit enables you to operate vehicle functions such as windscreen wipers, indicators, and headlights easily and safely, while at the same time keeping control of the steering wheel.

PARTNERS





FASTER BY DESIGN



SUNRISE MEDICAL







van**raam**® Let's all cycle









Innovating Vehicle Adaptations since 1983





otab

The leading car scheme for disabled people











Ellis PB is at New Brighton Sea Front.

22 Aug · Liverpool · 6

Thanks to Davinci Mobility for the chair that got Rachael around Italy, despite the Roman road builder's best efforts.

Gemma Partington I recommends Davinci Mobility.

Just back from a Saturday cycle along the Mersey coast (40k by myself sorry, Rob) - something I never imagined doing a year ago.

Cycling on my Davinci Mobility handcycle (nicknamed Vinny for Vincent Ross) over the last few months has been a godsend.

With the right infrastructure & support, disabled people can cycle.

I'm #PumpedUpToCycle, as Cycling UK like to say!



Lucy Chillery-Watson I recommends Davinci Mobility.

29 Aug 2019 · 🚱

If I could give more stars I would. Top customer service from beginning, middle and end. No job too big. My 6 year old daughter was involved in selecting her chair and specifications and they travelled over 4 hours to do a home assessment. Just amazing service. My daughter loves her Leggera Kids offroad wheelchair.



Cally Symms I recommends Davinci Mobility.

11 Apr 2019 · 🚱

a big thanks to neal and the team for helping me get my dream wheelchair

they were a great help from start to finish



Claire Watkins

12 May 2019 · 🕥

Davinci Mobility Neil Bevan thank you for arranging this, and making it possible! A LONG 26.2 miles complete in around ten hours from midnight to 10am throughout London. What a great company for sorting this out and donating their services for me to do this charity event. Amazing, thank you. If anyone would like to sponsor us please use the link below, if anyone would like information about mobility options available to make you independent please contact Davinci Mobility.



Kelly Rose I recommends Davinci Mobility.

27 Aug 2019 · 🕟

Thank you so much to Karl and Davinci for providing us with wheels...such a godsend and it has made our holiday. Thanks again for all your help



Simon Townley

10 Jul · 🕥

After 4 weeks in hospital I am FREE.... eg getting outside with the help of of my new trial rider from Davinci Mobility



Accesstheoutdoorguide is @ feeling thankful.

20 Aug 2019 · 🕥


On my way to the Edinburgh Fringe. I've been amazed how many people have stopped me to ask about my power wheel. My wheelchair and Trail Rider attachment is from Davinci Mobility. Taking the powerwheel attachment on the train has been stress free. Just detached it from the wheelchair and then it has tucked behind my chair for the journey. Edinburgh is so hilly that I need the extra power from the wheel to get about independently. Mark Townley




Bob Clauson I recommends Davinci Mobility.

I honestly can't rate DaVinci highly enough. I have known, and used, this company since it started out and they always provide a great service. They have just completed some modifications on my van, and knowing my circumstances and the need to have the van available, they fitted the jobs in around my free time, even collecting and delivering the van. I don't know if any other company who would be so obliging! I really appreciate the time and effort they put in to keeping me mobile. Whatever your mobility needs, DaVinci will have a solution.

The whole team are friendly, helpful and approachable. Thanks guys, you are all stars!













6-8 Carraway Rd, **Gilmoss Industrial Estate, Liverpool L11** 0EE

Opening times: 8:30 4:30 Monday to Friday



All Information and product is subject to change without notification. Please consult DaVinci Mobility for any queries you may have. If you are visually impaired, this document can be veiwed as a pdf at www.DaVinciMobility.co.





